

Reduce waste. Go green.

A reprocessing solution designed to:

- Meet sustainability goals
- Reduce hospital waste

HoverMatt® mattresses, Kendall SCD™ Sequential Compression Sleeves, Kendall DL™ Disposable Cable and Lead Wires and Nellcor™ pulse oximetry sensors are intended for single-patient use.

Fold and place HoverMatt® mattresses, compression sleeves, cables, lead wires and sensors in the provided collection hampers after the patient is discharged.

Kendall DL™ Disposable Cable and Lead Wires

After every patient use:



Detach all electrodes by pushing the pink button and releasing the head of the lead set from the electrode snap.



Disconnect the cable or lead set from the bedside monitor or telemetry box.



Leave the reusable adapter connected to the bedside monitor or telemetry box.

Kendall SCD™ Sequential Compression Sleeves

After every patient use:



Detach tubing from the sequential compression sleeves.



Leave tubing connected to the controller.



Finally, deposit used products in the hamper we have provided for reprocessing.



Do not use collection bins to discard the following items:



Sharps
No scalpels, needles, syringes, sutures or glass



Waste
No personal equipment, trash, patient care disposables or linens



Batteries
Remove from device but do not place in collection container



Kendall SCD™ controllers
Leave tubing connected to the Kendall SCD™ controllers



Kendall SCD™ tubing
Leave tubing connected to the Kendall SCD™ controllers



Reusable Kendall DL™ adapters
Leave adapters connected to bedside monitors or telemetry boxes

For more information call **844.268.2583**
or email **gogreen@cardinalhealth.com**

Sustainable Technologies™
A Cardinal Health business

© 2021 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, SUSTAINABLE TECHNOLOGIES, the Sustainable Technologies LOGO, KENDALL, KENDALL SCD and KENDALL DL are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. Lit. No. 2ST21-1649692 (09/2021)