

Reduce waste. Go green.

A sustainability solution designed to:

- **Meet recycling goals**
- **Reduce hospital waste**

Masimo™ pulse oximetry sensors, Kendall DL™ Disposable Cable and Lead Wires are intended for single-patient use.

Used Masimo™ pulse oximetry sensors, cables and lead wires should be placed in the provided collection bins after the patient is discharged.

After every patient use:
Kendall DL™ Disposable Cable and Lead Wires



Detach all electrodes by pushing the pink button and releasing the head of the lead set from the electrode snap.



Disconnect the cable or lead set from the bedside monitor or telemetry box.



Leave the reusable adapter connected to the bedside monitor or telemetry box.



Do not use collection bins to discard the following items:



Sharps
No scalpels, needles, syringes, sutures or glass



Waste
No personal equipment, trash, patient care disposables or linens



Batteries
Remove from device but do not place in collection container



Reusable Kendall DL™ adapters
Leave adapters connected to bedside monitors or telemetry boxes

For more information call **844.268.2583**
or email **gogreen@cardinalhealth.com**