

Device Collection Bins

Shipping Collection Bins to Provision

Small Collection Bin Preparation:

- To prevent damage to devices during shipping, pack the top of the small collection bin with packaging material.
- Small collection bins must be packaged inside a reinforced cardboard box and secured appropriately.
- Print a prepaid FedEx shipping label in Provision's SMART system or call Provision to request a label.

Standard Collection Bin Preparation:

- To protect the integrity of small, delicate devices, you can package such devices in other containers, such as specimen or denture cups, prior to placing them in the larger collection bin.
- To prevent damage to devices during shipping, pack the top of the collection bin with packaging material, such as bubble wrap, prior to closing the container lid.
- Close the lid of the collection container until you hear it snap in place and then lock the lid with the two zip ties which are taped to the underside of the lid.



Standard Collection Bin

24" x 16" x 9"

Small Collection Bin

12.5" x 7" x 4.5"

- Note that the standard collection bins are provided with shipping labels affixed to the container. When you are ready to ship the collection bin, you can either give the bin to any FedEx Ground driver on their normal delivery, or call Provision and we will arrange to have the container picked up.

DO NOT include any of the following items:

- Gowns
- Drapes
- Syringes
- Biopsy Forceps
- Devices with Pressure Cartridges
- Devices that contact Spinal Fluid
- Surgical Gloves
- Batteries
- Suture

340 SW Columbia Street, Bend, OR 97702
www.provisiondrs.com

541-508-7520