



Latex Classic

Delivers protection, dexterity and an advanced grip*

- Made with natural rubber latex
- 50 years of self-manufacturing expertise
- Suitable protection and dependable performance across a broad range of tasks
- Ideal outer glove when double-gloving or suitable as a standalone glove
- Interlocking, beaded cuff design reduces roll-down
- Anatomical fit and natural movement due to proprietary hand mold with an independent thumb design

The Protexis™ Latex Classic Surgical Glove delivers exceptional protection, dexterity and an advanced grip. Meets all relevant FDA and ASTM standards, including those for physical dimensions,[†] physical properties[‡] and freedom from holes.[‡] Documentation available upon request.

Clinical Application Matrix

Department	General	Cardiovascular	Dental/Maxillofacial	Ear, Nose and Throat (ENT)	Endovascular	Labor & Delivery	Laparoscopic/Robotics	Neuro	Obstetrics	Ophthalmology	Orthopedics	Pediatrics	Plastics	Thoracic	Urology	Vascular
Protexis™ Latex Classic Surgical Glove 2D72N55X-90X	●	●	●	●	●	●	●	●	●		●	●	●	●	●	●

This table was developed by a group of clinicians. It reflects current best practices of surgical glove usage per application. Ultimately, it is up to the discretion of the clinician to choose the right glove for the procedure.



Product information

			Thickness [†]								
Cat. no.	Size	Length	Finger	Palm	Cuff	Material	Color	Cuff type	Protein content [‡]	Qty/bx	Qty/cs
2D72N55X	5.5	11.5 in./ 292 mm	9.8 mil/ 0.25 mm	7.9 mil/ 0.20 mm	7.5 mil/ 0.19 mm	Natural rubber latex with nitrile polymer coating	 Cream	Beaded/ rolled	Less than 50 µg/ dm [‡]	50	200
2D72N60X	6										
2D72N65X	6.5										
2D72N70X	7	12 in./ 305 mm									
2D72N75X	7.5										
2D72N80X	8										
2D72N85X	8.5										
2D72N90X	9										

Properties (before aging)

Tensile strength (min)	≥ 24 MPa [†]
Stress at 500% elongation (modulus) (max)	≤ 5.5 MPa [†]
Ultimate elongation (elasticity) (min)	≥ 750% [†]
Puncture resistance (cuff) [§]	5N
Freedom from holes [‡]	0.65 AQL [‡]
Sterilization	Radiation
Accelerant	ZDBC (Zinc Dibutylthiocarbamate)



Protect what matters with Cardinal Health™ Protexis™ Surgical Gloves

Trust in high-quality gloves designed to protect both clinician and patient. Manufactured with strict quality controls and robust testing, you can depend on Protexis™ Surgical Gloves every time you enter the OR. Choose from a variety of comfortable, tactile options that can help support positive clinical outcomes. Rely on Protexis™ Surgical Gloves to help protect what matters most: you and your patients.

[†]As compared to Protexis™ Latex Surgical Gloves

[‡]In accordance with ASTM D 3577

[‡]Tested in accordance with ASTM D 5151

[‡]Protein content tested using ASTM D 5712

[§]Tested in accordance with AS/NZS 4179, min 5 N

Caution: This product contains natural rubber latex which may cause allergic reactions. Safe use of these gloves by latex-sensitized individuals has not been established.

© 2024 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO and PROTEXIS are trademarks of Cardinal Health and may be registered in the US and/or in other countries. Patent cardinalhealth.com/patents. Lit. No. 2MS19-1042127-02 (08/2024)

