



Product information

Cat. no.	Size	Length	Thickness [†]			Material	Color	Cuff type	Protein content [‡]	Qty/bx	Qty/cs
			Finger	Palm	Cuff						
2D72LS55	5.5	11.5 in./ 292 mm	9.8 mil/ 0.29 mm	9.1 mil/ 0.14 mm	7.1 mil/ 0.14 mm	Natural rubber latex coated with acrylic hydrogel polymer coating	 Translucent Yellow	Beaded/rolled	Less than 50 µg/dm [‡]	50	200
2D72LS60	6										
2D72LS65	6.5										
2D72LS70	7	12 in./ 305 mm									
2D72LS75	7.5										
2D72LS80	8										
2D72LS85	8.5										
2D72LS90	9										

Properties (before aging)

Tensile strength (min)	≥ 24 MPa [†]
Stress at 500% elongation (modulus) (max)	≤ 5.5 MPa [†]
Ultimate elongation (elasticity) (min)	≥ 750% [†]
Puncture resistance (cuff) [§]	5N
Freedom from holes [‡]	0.65 AQL [‡]
Sterilization	Radiation
Accelerant	ZDBC (Zinc Dibutylthiocarbamate)



Protect what matters with Cardinal Health™ Protexis™ Surgical Gloves

Trust in high-quality gloves designed to protect both clinician and patient. Manufactured with strict quality controls and robust testing, you can depend on Protexis™ Surgical Gloves every time you enter the OR. Choose from a variety of comfortable, tactile options that can help support positive clinical outcomes. Rely on Protexis™ Surgical Gloves to help protect what matters most: you and your patients.

[†]As compared to Protexis™ Latex Surgical Gloves

[‡]In accordance with ASTM D 3577

[‡]Tested in accordance with ASTM D 5151

[‡]Protein content tested using ASTM D 5712.

[§]Tested in accordance with AS/NZS 4179, min 5 N

Caution: This product contains natural rubber latex which may cause allergic reactions. Safe use of these gloves by latex-sensitized individuals has not been established.

© 2024 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO and PROTEXIS are trademarks of Cardinal Health and may be registered in the US and/or in other countries. Patent cardinalhealth.com/patents. Lit. No. 2MS19-1042514-03 (08/2024)

